


# microVISC PRO

Easy to Use, Portable, Small Sample Viscometer



## Key Benefits



Easy to use measurements with built-in touchscreen controls.



Rapid measurements in less than 1 minute.



Measure Newtonian and Non-Newtonian samples- 100µL



Test at multiple temperatures with fast built-in temperature control.

The portable microVISC PRO is your solution for quick accurate viscosity measurements anywhere. Any user can run samples based on previously set parameters.

microVISC PRO employs RheoSense's patented VROC® technology and offers repeatable viscosity measurements. Using samples as small as 100 microliters, VROC saves time and resources.

## Specifications:

- Minimum Sample Volume –100 µL
- Shear Rate Range, s<sup>-1</sup>: 1 ~ 8,300
- Viscosity Range: 0.3 ~ 650,000 mPas
- Temperature range 15 °C – 60 °C
- Accuracy up to 2% of Reading
- Repeatable Measurement up to 0.5% of Reading
- 21 CFR Part 11 Compliant
- Additional configurations available